

easy sterile connections for high flow applications, even in non-sterile environments. AseptiQuik X's "TWIST-PULL-TWIST" design enables users to quickly transfer large volumes of media easily with less risk of operator error. The connector's robust design provides reliable performance without the need for clamps, fixtures or tube welders. Biopharmaceutical manufacturers can now make 1" sterile connections with the quality and market availability they expect from the leader in single-use connection technology.

FEATURES	BENEFITS
TWIST-PULL-TWIST Design	Intuitive three-step connection process reduces risk of operator error
Membrane pull tabs	Ensure simultaneous and secure removal of both membranes
Robust construction	Repeatable and reliable performance with no additional hardware required
Integrated lock ring	Secures final connection preventing disassembly
CPC Click	Audible confirmation of completed assembly steps

Specifications • • •



PRESSURE:

Up to 60 psi, 4.1 bar

TEMPERATURE:

39°F to 104°F (4°C to 40°C)

STERILIZATION:

Gamma: Up to 50kGy irradiation

Autoclave High Temp (HT) Version: Up to 266°F (130°C) for 30 minutes

TERMINATION SIZES:

3/4" and 1" ID hose barb (19.0mm and 25.4mm). 1-1/2" sanitary

MATERIALS:

Main Components:

PVDF (white), USP Class VI, ADCF

Lock Ring:

PVDF (blue), USP Class VI, ADCF

Polycarbonate (blue, standard version), USP Class VI, ADCF

Polycarbonate (white, HT version), USP Class VI, **ADCF**

Polypropylene (clear), USP Class VI, ADCF


Seals:

Silicone (clear), platinum-cured, USP Class VI, ADCF

Membrane:

Polyethylene (standard version), USP Class VI, ADCF Hydrophobic polyethersulfone (HT version), USP Class VI, ADCF, PTFE strip sticker

This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area.





ASEPTIQUIK® X SERIES DIMENSIONS



B = Total Length

= Hose Barb Length





Coupling Bodies

POLYCARBONATE with Blue Pull Tabs



TERMINATION 3/4" HOSE BARB 1" HOSE BARB 1-1/2" SANITARY PART NO. AQX17012 AQX17016 AQX33024 3.54 3.54 3.54

3.93 4.04 4.18 1.50 1.53 1.73

POLYCARBONATE HT with White Pull Tabs



TERMINATION 3/4" HOSE BARB 1" HOSE BARB 1-1/2" SANITARY PART NO. AQX17012HT AQX17016HT AQX33024HT 3.54 3.54 3.54 3.93 4.04 4.18 1.50 1.53

Coupling Inserts

• POLYCARBONATE with Blue Pull Tabs



TERMINATION 3/4" HOSE BARB 1" HOSE BARB 1-1/2" SANITARY

PART NO.AQX22012
AQX22016
AQX44024

A 3.76 3.76 3.76 4.34 4.45 1.50 1.53


POLYCARBONATE HT with White Pull Tabs



TERMINATION 3/4" HOSE BARB 1" HOSE BARB 1-1/2" SANITARY PART NO. AQX22012HT AQX22016HT AQX44024HT 3.76 3.76 3.76

4.34 4.45 4.58 1.50 1.53 1.74

TWIST-PULL-TWIST Assembly Procedure



Align male and female connectors using the START alignment feature of the blue lock ring with the rib indicator of the white body. Insert the two halves together.



Twist the blue lock ring clockwise (approximately 75°) until audible "CPC Click" is heard. Alignment of the small blue lock ring rib indicator with the white body's rib indicator confirms final rotation is complete.



Snap the membrane pull tabs together and pull from connector.



Twist the blue lock ring clockwise until the final audible "CPC Click" is heard (90° from the initial starting point). Alignment of the long blue lock ring rib indicator with the white body's rib indicator confirms initial rotation is complete.